

Ural Instrument
Engineering Plant JSC



МОБИВЕНТ ОКСИ



+7 (343) 359-93-76

mail@upz.ru

www.upz.ru



Designed for high-flow oxygen therapy in adult and pediatric patients, with a built-in humidification system

Monitoring of temperature and humidity of the breathing mixture, flow and pressure in the patient's circuit, blood oxygen saturation and heart rate, graphic monitoring;

Modes of non-invasive ventilatory support through a mask or nasal cannula;

Extended alarm system for ventilation parameters and technical alarms;

Informative color 3.5-inch touch-sensitive display;

Ultrasonic flow and oxygen sensor;

SpO2 module.

